



ANALYTICAL PROCEDURES

For DR/2000 and DR/3000 Instruments


Method 10023

NITROGEN, AMMONIA, Low Range, Test 'N Tube™ *

(0 to 2.50 mg/L NH₃-N)

Salicylate Method**

For water, wastewater, and seawater



1. Enter the stored program for Low Range Test 'N Tube Nitrogen, Ammonia.

Press: **3 4 2 READ/ENTER**


The display will show:

Dial nm to 655

Note: See Instrument Setup on page 3 to enter this method into the DR/2000.

Note: DR/2000's with software versions 3.0 and greater will display **P** and the program number.


Note: DR/2000's with software versions 3.0 and greater will not display **DIAL TO** message if the wavelength is already set correctly. The display will show the message in step 3. Proceed with step 4.



2. Rotate the wavelength dial until the display shows:

655 nm


Note: For DR/3000 instruments, set the wavelength to 655 nm and press **CLEAR**.



3. Press: **READ/ENTER**

The display will show:

mg/L N Vial LR




4. Remove the caps from 2 AmVer™ Diluent Reagent LR vials. Add 2 mL of deionized water to 1 vial (the blank). Add 2 mL of sample to the other (the sample).

Note: If sample cannot be analyzed immediately, see Sampling and Storage on page 7.


* AmVer™, PourRite™, TenSette®, Test 'N Tube™, and Voluette™ are trademarks of Hach Company.

** Adapted from Clin. Chim. Acta 14: 403 (1966).


NITROGEN, AMMONIA, Low Range, Test 'N Tube™, continued



5. Using a funnel, add the contents of 1 Ammonia Salicylate Reagent Powder Pillow for 5 mL Sample to each vial.




6. Using a funnel, add the contents of 1 Ammonia Cyanurate Reagent Powder Pillow for 5 mL Sample to each vial.



7. Cap the vials tightly and shake thoroughly to dissolve the powder.


Note: A green color will develop if ammonia is present.



8. Press: **SHIFT TIMER**

A 20-minute reaction period will begin.

Note: For DR/3000's, press **20 TIMER**.




9. When the timer beeps the display will show:

mg/l N Vial LR


Place the COD Vial Adapter into the cell holder with the marker to the right.

Note: For DR/3000s, the groove in the adapter faces the front of the instrument.




10. Clean the outside of both vials with a towel.

Note: Wipe with a damp cloth and follow by a dry one to remove fingerprints and other marks.



11. Place the blank into the vial adapter with the Hach logo facing the front of the instrument. Place the cover on the adapter.



12. Press: **ZERO**

The display will show:


WAIT

then:


0.00 mg/l N Vial LR

Note: For DR/3000's, press **MANUAL PROGRAM**. Then press **ZERO**. Press **ZERO** again if the display does not show **0.000**. Enter the concentration factor (**1.502**) and press **CONC FACTOR**. Then press **2 CONC**.

NITROGEN, AMMONIA, Low Range, Test 'N Tube™, continued



13. Place the prepared sample into the vial adapter with the Hach logo facing the front of the instrument. Place the cover on the adapter.



14. Press: **READ/ENTER**

The display will show:

WAIT

then the result in mg/L ammonia nitrogen (N) will be displayed.

Note: In the constant on mode, pressing **READ/ENTER** is not required. **WAIT** will not appear. When the display stabilizes, read the result.

Note: The result can be expressed as mg/L ammonia (NH₃) by multiplying the mg/L N result by 1.22.

Instrument Setup

DR/2000 with 1.261 or 1.27 software

Enter the calibration as an operator-programmed calibration. Follow the steps in the Instrument Operation section of the *Instrument Manual*. Store the method as follows:

nm = 655

Decimal = 00.00

Units = mg/l

Symbol = N Vial LR

Timer 1 = 20:00

Enter the calibration with 0.000 absorbance values for zero and #1 standard. To do this, do not place anything in the sample cell compartment. Begin by storing standards 0 and 1 as the concentrations shown in the table below. Accept 0.000 Abs as the absorbance value for all standards. Store the calibration values by pressing **SHIFT READ/ENTER**.

NITROGEN, AMMONIA, Low Range, Test 'N Tube'™, continued

Next, edit the absorbance values for the standards to the values given below.
Follow the steps given in the Instrument Operation section of the *Instrument Manual*.

Std	Conc	Abs
#0	0.00	0.000
#1	3.00	1.968

The method is now stored as an operator-programmed method number between 950 and 999. Record the method number for future reference when using this method.

DR/2000 with software version 2.0 or 2.2




Enter the calibration as an update to Hach-stored programs.

1. Press: 

1. Press:  

2. Press:    

3.

4. Within 3 seconds, press:   

The display will show: **ENTER nm**

5. Press:    

Note: If you make an error, press **SHIFT CLEAR** and re-enter the number. When the number is correct, press **READ/ENTER**.

The display will show: **DECIMAL? 00.00**

6. The decimal point is correctly positioned. Press **READ/ENTER**. The display will show: **UNITS?**

7. Use the arrow keys to select the appropriate unit of measure. Press the down arrow key twice. The display will show: **mg/l**

8. Press **READ/ENTER** when the correct unit of measure is displayed.
The display will show: **SYMBOL?**

9. Construct the correct symbol display: **N Vial LR**

- Select letters and regular numbers by scrolling to the correct symbol with the arrow keys.
- To make a letter uppercase, press the **SHIFT** key.
- The space is the character displayed after one press of the down arrow key.
- Accept each symbol by pressing **READ/ENTER**.

NITROGEN, AMMONIA, Low Range, Test 'N Tube™, continued

- e. To end symbol entry, press **READ/ENTER** a second time after accepting the last character.

When the instrument is out of symbol entry mode, the display will show: **TIMER?**

10. This method has 1 timed step, so press **SHIFT TIMER**. The display will show:
MM:SS TIME 1?

11. Enter a timer value of 20 minutes. Press:



12. Press **READ/ENTER** to accept the timer value. The display will show:
MM:SS TIME 2?

13. Press **READ/ENTER** to complete the timer entry. The display will show: **#1 Data**

14. Enter the following 12 numbers as shown. Complete each number entry by pressing the **READ/ENTER** key.

Display	Number Entry
#1 Data	0
#2 Data	4883
#3 Data	4883
#4 Data	4883
#5 Data	4883
#6 Data	4884
#7 Data	4883
#8 Data	4883
#9 Data	4883
#10 Data	10922
#11 Data	512
Checksum	15037

The final number is a check value that determines if the data sequence was entered correctly. If an error was made during number entry, the display will return to the prompt for #1 Data and the entire sequence must be re-entered. If all numbers are correctly entered, the display will return to the method prompt and is ready for use.

DR/2000 with software version 3.0 or 3.1

1. Turn the instrument on. Press **SHIFT METHOD** to enter the configuration mode. The display will show: **MOMENTARY** or **CONSTANT ON**
2. Press the up arrow key twice to select **HACH UPDATE**. Press **READ/ENTER**. The display will show: **ENTER #:**

3. Press:

The display will show: **P342 ENTER nm**

4. Press:






NITROGEN, AMMONIA, Low Range, Test 'N Tube'TM, continued

Note: If you make an error, press **SHIFT CLEAR** and re-enter the number. When the number is correct, press **READ/ENTER**.

The display will show: **P342 Decimal? 00.00**

5. The decimal is correctly positioned. Press: **READ/ENTER** to accept the correct position. The display will show: **P342 UNITS?**
6. Use the arrow keys to select the appropriate unit of measure. Press the down arrow key twice. The display will show: **P342 mg/l**
7. Press **READ/ENTER** when the correct unit of measure is displayed. The display will show: **P342 mg/l**
8. Construct the display to read the correct symbol. The symbol must be entered exactly as shown including dashes and spaces between characters: **N Vial LR**
 - a. Select letters and numbers by scrolling to the correct character with the arrow keys.
 - b. To make a letter uppercase, press the **SHIFT** key.
 - c. The space is the character displayed after one press of the down arrow.
 - d. Make sure to enter the display line EXACTLY as shown, including the spaces. Do not enter trailing spaces.
 - e. Accept each symbol by pressing **READ/ENTER**.
 - f. When the last character of the symbol is accepted with the **READ/ENTER** key, press **READ/ENTER** a second time to end display entry mode.

When the instrument is out of symbol entry mode, the display shows: **P342 TIMER?**

9. This method has 1 timed step, so press **SHIFT TIMER**. The display will show: **MM:SS TIME 1?**
10. Enter a timer value of 20 minutes. Press:     
11. Press **READ/ENTER** to accept the timer value. The display will show: **MM:SS TIME 2?**
12. Press **READ/ENTER** to complete the timer entry. The display will show: **# 0 STANDARD**
13. Press **READ/ENTER** to display the zero data pair. The display will show: **0.000 Abs 00.00 mg/l**
14. Press **READ/ENTER**. The display will show: **#1 STANDARD**
15. Press **READ/ENTER**. The display will prompt for entry of the first concentration point: **#1 00.00 mg/l**
16. Enter concentration point #1 from the table below by pressing **0300** so that the display shows: **# 1 03.00 mg/l**

NITROGEN, AMMONIA, Low Range, Test 'N Tube™, continued

17. Press **READ/ENTER**. The display will prompt for entry of the first absorbance point:
1 0.000 Abs

18. Enter the absorbance point #1 from the table below by pressing **1968** so that the display shows: **# 1 1.968 Abs**

19. Press **READ/ENTER**. The display will show the first data pair: **1.968 Abs 03.00 mg/l**

20. Press **READ/ENTER** to accept the first data pair. The display will show:
#2 STANDARD

The data pair values from the table below are now entered.

Standard	Concentration	Absorbance
#0	[0.00] mg/l	[0.000] Abs
#1	[3.00] mg/l	[1.968] Abs

When the last data pair is entered the display will show: **#2 STANDARD**

21. Press **SHIFT READ/ENTER** to complete data point entry. The display will show: **#:**

22. Enter the validation number: **5164** so that the display shows: **#: 5164**

23. Press **READ/ENTER**. The display will show: **COMPLETED**

then: **P342 mg/l N Vial LR**

Note: If the display shows **INCORRECT #**, then prompts again for the validation number, you may have made an error during data entry. Make sure the validation number is correct. If so, then the error occurred during some other portion of the method entry. Press **METH** and respond to the **ABORT?** message by pressing **READ/ENTER**, then re-enter the method.

The instrument is now ready for use with Method 342.

Sampling and Storage

Collect samples in clean plastic or glass bottles. Best results are obtained with immediate analysis. If chlorine is known to be present, add 1 drop of 0.1 N Sodium Thiosulfate Standard Solution for each 0.3 mg/L Cl_2 in a 1 liter sample. Preserve the sample by reducing the pH to 2 or less with concentrated Hydrochloric Acid (at least 2 mL). Store at 4 °C (39 °F) or less. Preserved samples may be stored up to 28 days. Before analysis, warm samples to room temperature and neutralize with 5.0 N Sodium Hydroxide Standard Solution. Correct the test result for volume additions.

Accuracy Check

Standard Additions Method

- Snap the neck off a Nitrogen, Ammonia Standard Solution Ampule, 50 mg/L $\text{NH}_3\text{-N}$.
- Use the TenSette® Pipet to add 0.1, 0.2, and 0.3 mL of standard to three 25 mL samples. Mix thoroughly.

NITROGEN, AMMONIA, Low Range, Test 'N Tube'™, continued

- c. Analyze each sample as described above. The nitrogen concentration should increase 0.20 mg/L for each 0.1 mL of standard added.
- d. If these increases do not occur, see Standard Additions in the *Procedures Manual* for more information.

Standard Solution Method

To check accuracy, use the Nitrogen, Ammonia Standard Solution, 1.0 NH₃-N mg/L listed under Optional Reagents. Or, dilute 1 mL of solution from a 50 mg/L Voluette™ Ampule Standard for Nitrogen, Ammonia to 50 mL with deionized water using a 50 mL volumetric flask.

Precision

DR/2000: In a single laboratory, using a standard solution of 1.5 mg/L ammonia nitrogen (NH₃-N) and 2 representative lots of reagent with the DR/2000 Spectrophotometer, a single operator obtained a standard deviation of ± 0.03 mg/L N.

DR/3000: In a single laboratory, using a standard solution of 1.5 mg/L ammonia nitrogen (NH₃-N) and 2 representative lots of reagent with the DR/3000 Spectrophotometer, a single operator obtained a standard deviation of ± 0.03 mg/L N.

Interferences

The following ions may interfere when present in concentrations exceeding those listed below:

Interfering Substance	Interference Level and Treatment
Calcium	2500 mg/L as CaCO ₃
Iron	Determine the amount of iron present in the sample following one of the total iron procedures. Add the same iron concentration to the deionized water in <i>step 4</i> . The interference will then be successfully blanked out.
Magnesium	15,000 mg/L as CaCO ₃
Nitrite	30 mg/L as NO ₂ ⁻ -N
Nitrate	250 mg/L as NO ₃ ⁻ -N
Orthophosphate	250 mg/L as PO ₄ ³⁻ -P
pH	Acidic or basic samples should be adjusted to about pH 7. Use 1.0 N Sodium Hydroxide Standard Solution for acidic samples and 1.0 N Hydrochloric Acid Standard Solution for basic samples.
Sulfate	300 mg/L as SO ₄ ²⁻
Sulfide	Will intensify the color. Measure about 350 mL of sample in a 500 mL Erlenmeyer flask. Add the contents of one Sulfide Inhibitor Reagent Powder Pillow. Swirl to mix. Filter the sample through a folded filter paper. Use the filtered solution in <i>step 4</i> .
Other	Less common interferences such as hydrazine and glycine will cause intensified colors in the prepared sample. Turbidity and color will give erroneous high values. Samples with severe interferences require distillation. Hach recommends the distillation procedure using the Hach General Purpose Distillation Set. See <i>OPTIONAL APPARATUS</i> at the end of this procedure.

NITROGEN, AMMONIA, Low Range, Test 'N TubeTM, continued

Summary of Method

Ammonia compounds combine with chlorine to form monochloramine. Monochloramine reacts with salicylate to form 5-aminosalicylate. The 5-aminosalicylate is oxidized in the presence of a sodium nitroprusside catalyst to form a blue-colored compound. The blue color is masked by the yellow color from the excess reagent present to give a green-colored solution.

NITROGEN, AMMONIA, Low Range, Test 'N TubeTM, continued

Reagents and Apparatus

REQUIRED REAGENTS

AmVer TM Reagent Set for Nitrogen, Ammonia, Low Range (50 vials)	26045-45
Includes: (1) 23952-66, (1) 23954-66, (1) 272-42, (50) AmVer TM Diluent LR Vials*	

Description	Unit	Cat. No.
AmVer TM Diluent Reagent, Test 'N Tube TM low range vials	50/pkg	*
Salicylate Reagent Powder Pillows, 5 mL Sample	50/pkg	23952-66
Cyanurate Reagent Powder Pillows, 5 mL Sample	50/pkg	23954-66
Water, deionized.....	100 mL	272-42

REQUIRED APPARATUS

Description	Quantity Required		Cat. No.
	Per Test	Unit	
COD Vial Adapter, DR/2000 and DR/3000	1	each	44799-00
Funnel, micro (for reagent addition).....	1	each	25843-35
Pipet, TenSette [®] , 0–10 mL.....	1	each	19700-10
Pipet Tips for 19700-10	1	50/pkg	21997-96
Test Tube Rack.....	1–3	each	18641-00
Safety Bulb	1	each	14651-00

OPTIONAL REAGENTS

Hydrochloric Acid, ACS.....	500 mL	134-49
Nitrogen, Ammonia Standard Solution, 1.0 mg/L NH ₃ –N.....	500 mL	1891-49
Nitrogen, Ammonia Standard Solution, 50 mg/L NH ₃ –N, 10 mL ampules	16/pkg	14791-10
Nitrogen, Ammonia Standard Solution, 2 mL PourRite TM ampules, 50 mg/L NH ₃ –N	20/pkg	14791-20
Sodium Hydroxide Standard Solution, 5.0 N	50 mL SCDB	2450-26
Sodium Hydroxide Standard Solution, 1.0 N	100 mL MDB	1045-32
Sodium Thiosulfate Standard Solution, 0.1 N	100 mL SCDB	323-32
Sulfide Inhibitor Powder Pillows	100/pkg	2418-99
Sulfuric Acid, 1.0 N.....	10 mL MDB	1270-32
Water, deionized.....	4 L	272-56

OPTIONAL APPARATUS

Ampule Breaker Kit.....	each	21968-00
Ampule Breaker Kit, PourRite TM	each	24846-00
Distillation Apparatus Set, general purpose	each	22653-00
Filter Paper, folded.....	100/pkg	1894-57
Flask, Erlenmeyer, 500 mL.....	each	1082-49
Flask, volumetric, 50 mL	each	14574-41
Funnel, analytical (for filtering).....	each	1083-68
Heater and Support Apparatus (for distillation), 115 Vac	each	22744-00
Heater and Support Apparatus (for distillation), 230 Vac	each	22744-02
pH Indicator Paper, 1 to 11 pH	5 rolls/pkg	391-33
Thermometer, -20 to 110 °C	each	566-01
Thermometer, -10 to 260 °C	each	20959-26

* This item is not sold separately. Please order the complete set (cat. no. 26045-45) as a replacement.



HACH COMPANY
WORLD HEADQUARTERS
P.O. Box 389
Loveland, Colorado 80539-0389
Telephone: (970) 669-3050
FAX: (970) 669-2932

HACH EUROPE
Chaussée de Namur, 1
B-5150 Floriffoux (Namur), Belgium
Telephone: (32)(81) 44.71.71
FAX: (32)(81) 44.13.00

FOR TECHNICAL ASSISTANCE, PRICE INFORMATION AND ORDERING:

In the U.S.A. - **Call toll-free 800-227-4224**

Outside the U.S.A. - **Contact the HACH office or distributor serving you.**

On the Worldwide Web - **www.hach.com**; E-mail - **techhelp@hach.com**
