Sonda 50 48

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
50 48 probe
Refer to the meter manual for more information.

Caution

⚠ Personal injury hazard. Glass components can break. Handle with care to prevent cuts.

⚠ Chemical exposure hazard. Refer to the current material safety data sheets (MSDS) for safety protocols.

Introduction

This probe is a combination gel-filled electrode with a built-in temperature sensor for pH, ORP, and conductivity measurements in water treatment, galvanic applications and pool maintenance.

Technical information

Range	0 to 14 pH; ± 2000 mV; 20 µS/cm to 300 mS/cm
Junction	Ceramic
Reference	Ag/AgCl
Temperature	Pt 1000
Electrodes	EC, ORP: platinum
Minimum depth	10 mm (0.4 in.)
Operating temp	0-80 °C (32-176 °F)
Limits	Colloids, HF
Connector type	MP-8

Preparation

First time use: put the probe in ethanol for 15 s. Rinse with deionized water.

Connect the probe

⚠ Do not turn or twist the probe connector.

1 Calibration

Caution

⚠ Air bubbles under the probe tip when submerged can cause slow stabilization or error in measurement. If bubbles are present, gently shake the probe until bubbles are removed.

For 1 or 2-point pH calibration: push ✓ after 1 or 2 standards.

For conductivity, use 1-point calibration and select the standard nearest to the expected measurement range.

2 Measurement

Continuous: push ✓ 2 times.

Change the parameter: push and hold ✓ during stabilization.

Maintenance and storage

Cleaning: soak the probe in a cleaning solution. Rinse with deionized water.

Hydration: soak the dry probe in a dilute HCl solution for several hours.


Storage: put the probe in the protector with conductivity standard 1413 µS/cm.

Troubleshooting

Be sure that the correct probe type is used. Do not hold the tubes by hand to prevent temperature changes.

Warranty

6 months for manufacturing faults only. Damage from use is not covered.

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