

Item #	Product Description - Check List	Qty/Kit
2943200	Xenosep Starter Kit with Solvent Recovery	1
	Each Kit Contains the Following Items:	
2943222	Tube, Collection, 30x200mm	1
2943224	Clamp, Plastic, TS 24/25	3
2943213	Flask, Flat Sided Erlenmeyer, TS 24/25, 125mL	7
2943212	Flask, Round Flat Bottom, TS 24/25, 125mL	1
2943215	Condenser, Crossover w/#2HB's	1
2943214	Eluter	1
2943225	Stopcock, Glass, 120° 3 way, 4mm ID	1
2943701	Standard, EPA 1664A, 1 x 30mL, 4mg/mL	5
2943229	O-Ring, Collection Tube	1
	Clamp, Tubing, Acetal	1
	Stand, Tube	1
2943217	Funnel, 300 mL, 47mm	1
2943216	Holder, Starter, 47mm	1
2943223	O-Ring, Coupler, 47mm	2
2943218	Coupler, PTFE, 47mm	1
2943219	Support, Stainless Steel, 47mm	1
2943220	Clamp, Aluminum, 47mm, Wide Opening	1
2943211	Flask, TS40/35, 1L	1
2943221	Kit, Tubing, Tygon(R), 3/8"IDx5/8"OD, 5'	1
2943301	Filters, SPE, 1664A, 47mm, XM	24
2943401	Pipets, Transfer, Polyethylene, 15mL, XM	24
2943501	Columns, Sodium Sulfate, 10g, XM	24
2943237	Pans, Aluminum, 105mm, XM	24
	Method Summary	1
	Operating Instructions	1
	CD - Video Demonstration	1